



'Expert mode'



Vladimir Mikryukov

UT International Summer University, Tartu, August 01-04 2022



https://mycology-microbiology-center.github.io/Metabarcoding2022/

https://t.ly/dvZS

Metabarcoding: from Lab to Bioinformatics

Metabarcoding: from Lab to Bioinformatics (UT International Summer University, 2022)

Metabarcoding: from Lab to Bioinformatics

University of Tartu, 2022


Data used during the course

"Expert mode" commands


naividual projects


Slides (will be released after the course)

Expert mode in PipeCraft2



Docker






Windows Terminal

+ Windows Subsystem for Linux (WSL)





https://docs.microsoft.com/en-us/windows/wsl/install


```
cd ~ cd /data
```

```
cd /data/OC
cat QC.md5
md5sum -c QC.md5
md5sum *.fq.gz > Checksum.md5
md5sum MiSeq R\{1,2\}.fq.qz
```

```
less QC.md5
```

```
grep "MiSeq" QC.md5
grep "MiSeq" QC.md5 | wc -1
```

```
zgrep "GTGAATCATCGAATCTTTG" \
   --color=always \
   "MiSeq R1.fq.gz" \
   l less -R
zgrep "GTGAATCATCGAATCTTTG" \
   "MiSeq R1.fq.gz" \
   | grep -v "TAAGCGGCGGACT" \
   sort | uniq \
   > some sequences.txt
```





USEARCH


https://drive5.com/usearch/



Robert Edgar







Metabarcoding: from Lab to Bioinformatics

Metabarcoding: from Lab to Bioinformatics (UT International Summer University, 2022)

Metabarcoding: from Lab to Bioinformatics

University of Tartu, 2022

Data used during the course

"Expert mode" commands

naividual projects

Slides (will be released after the course)

Docker images

Working docker images

Image	Software
ewels/multiqc:latest	mutliqc v1.12
staphb/fastqc:0.11.9	fastqc v0.11.9
pipecraft/cutadapt:3.5	cutadapt v3.5, seqkit v2.0.0, python3, biopython
pipecraft/dada2:1.20	dada2 v1.20, seqkit v2.0.0, lulu v0.1.0, R
pipecraft/reorient:1	fqgrep v0.4.4, seqkit v2.0.0, mothur v.1.43.0
pipecraft/trimmomatic:0.39	trimmomatic 0.39, seqkit v2.0.0
pipecraft/vsearch:2.18	vsearch v2.18, seqkit v2.0.0, GNU parallel
pipecraft/itsx:1.1.3	ITSx v1.1.3, seqkit v2.0.0, mothur v1.46.1
pipecraft/blast:2.12	BLAST v2.12.0+
pipecraft/deicode:0.2.4	DEICODE v0.2.4, qiime2-2002.2
pipecraft/fastp:0.23.2	fastp v0.23.2
pipecraft/blast:2.12	BLAST 2.12.0+, biopython, python3, gawk

https://pipecraft2-manual.readthedocs.io/en/latest/user_guide.html#working-docker-images

Run VSEARCH container

```
docker run -it
  -v /c/Metabarcoding:/data
  pipecraft/vsearch:2.18
```

I. Dereplicate sequences

```
cd /data
```

```
vsearch
```

- --derep fulllength MC 1 .fasta
- --output MC 1 dereplicated.fasta
- --fasta width 0
- --relabel_sha1
- --sizeout

II. UNOISE clustering

```
vsearch
```

- --cluster_unoise MC_1_dereplicated.fasta
- --strand both
- --minsize 4
- --unoise alpha 2
- --fasta width 0
- --sizein --sizeout
- --centroids MC_1_unoise.fasta

Count sequence numbers

```
grep -c "^>" MC_1_dereplicated.fasta
grep -c "^>" MC 1 unoise.fasta
```

Loop through all samples

```
for i in
chimera_Filtered_out/*.fasta; do
  echo ${i}
```

done

Loop through all samples

```
vsearch
  --derep_fulllength ${i}
  --output -
| vsearch
  --cluster_unoise -
  --centroids $(basename ${i})
```